


Dear Amira Customer,

You can install the software from the DVD included in this package. Once successfully installed, Amira has to be activated with a license key. In order to issue a license key to you, we need the Computer ID of your workstation and the License Serial Number of each product you desire to activate. Each product is locked to a specific computer ID, this means that this product will execute only on the system(s) for which it is licensed.

How can I find the Computer ID?

In order to obtain the computer ID please start Amira and choose Help/License Manager. If no license is installed, the license manager pops up automatically upon the start. The License Manager will show your Computer ID in the rightmost text field (See image below, field "Computer ID").



Where can I find the License Serial Number?

Each software license you purchased has a License Serial Number indicated in the Delivery Note you received via email. Please check the Delivery note with attention, you will find a serial number starting with 257- for each license option. This serial number is the License Serial Number.

What should I do to permanently activate Amira?

Please go to <http://www.amira.com/lkr> and provide this Computer ID along with your License Serial Number(s) and contact information.

If you have any problems, please contact us at

North America, Central & South America

Visage Imaging, Inc.
1815 Aston Avenue - Suite 107
Carlsbad, CA 92008-7310, USA
Phone toll-free (US only): +1-866-627-1673
Phone: +1-760-494-9698
amirasupport@visageimaging.com

Europe, Africa, Middle East & Asia / Pacific

Visage Imaging GmbH
Lepsiusstr. 70
12163 Berlin, Germany
Phone: +49 30 700 968 - 20
amirasupport@visageimaging.com